**MULTI  
CAP**

**FLAT  
& DOC**

**BASE**

**SOFT**

**CUP**

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Kopfbänder

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Quick Insert Elektrode

This user manual provides the correct usage and maintenance of the MultiCap EEG system-caps in the versions "Flat", "Base", "Soft" and "Cup", with flexible positionable electrodes and their accessories.

## Legend of symbols

Symbol	Meaning
	Attention! read this instruction aware
	Serial number
	Not dry-cleaning!
	Washing careful, minimum rpm!
	Not dashing!
	Not ironing!

## Warnings and indications

- English**
-  By use of this product, keep this indications necessarily.
  -  By improperly use, the guarantee can get lost.
  -  During the recording, mobile phones, radios, other transmitters or receivers, x-ray devices, HF-devices and/or strong electrical appliances should be avoided.
  -  Electrodes of the MultiCap- Flat, MultiCap-Base and the MultiCap-Soft do not use nearly nuclear magnetic resonance imaging (NMRI).
  -  Direct sunlight should be avoided.
  -  Excessive stretching, twisting or tangling of the lead wires should be avoided. The cables should be stored hanging.
  -  By reuse with several patients, the product must be cleaned and disinfected between the recordings, in order to avoid infection-transfers (see page 15).
  -  Use only special pastes and gels during the recording, which are recommended by the manufacturer (see page 13).
  -  Check the functional ability of the product and the accessories sequentially in regular distances to guarantee an undisturbed application.
  -  Electrodes transfers data at a connected EEG device. If the product should take a damage through a incorrect data transport by the hard- or software, the user or the manufacturers of the hard- or software are be liable.
  -  Electronical sensors or electrodes are breakable and must be treated carefully. For this, lay down the electrodes for example on a towel anytime. Do not touch metal or metallic objects with the electrode

## Designated use

MultiCaps are textile EEG recording caps with flexible positionable electrodes, which are pre-positioned by the international 10-20 recording system. With MultiCaps it can be used different types of electrodes for the recording. MultiCaps can be used in clinical practices, with neurological or functional diagnostics and also for scientific researches.

MultiCaps are determined for use with Encephalographs together. The system caps are made from elastic textile material. That provides a very good positioning of all electrodes without necessary additional adaptation. Accounting of the flexibility of the material, the required number of EEG system-caps will be reduced to a minimum, for the most different head sizes.

You can get the MultiCap-Base as single system cap or as Full-System. The full-system contains also required products for the EEG recording.

## Specification

Parameter	Wert
Electrode connectors	1,5mm DIN ( <i>MultiCap-Base and Multi Cap-Cup</i> ), 25pole multi connector with adapter ( <i>MultiCap-Soft</i> )
Length of electrode wire	120cm +/- 5cm ( <i>MultiCap-Base</i> ), 150cm +/- 5cm ( <i>MultiCap-Soft and MultiCap-Cup</i> )
Electrode polarization	<= 100 mV
Electrode impedance	< 5 kOhm
Weight of electrode	5g +/- 0,5g
Weight of cap without electrodes	Max 100g
Mean time between failures (MTBF)	Up to 500 cycles* of use
Working temperature	+10°C bis +45°C
Storage temperature	+5°C bis +45°C
Transportation temperature	1-50°C bis +50°C
Humidity by storage	Up to 97% without condensation

\* Cycles refer of recordings, included the time for washing and disinfection of the cap and the electrodes.

## Procedure

1. Unpack the product and check the completeness and appearance of all parts. In case of missing or damaged parts please contact your supplier.
2. Fill the syringe with gel and install the supplied blunt needle. The needle should sit tightly on the syringe.

English

### Attention!

Use only the supplied blunt needle of plastic. A sharp needle of metallic can damage the skin of the patient.

## Instruction for use

1. Measure the head circumference of your patient and choose the proper cap size.
2. Pull the cap on the head of the patient. Beginning from the forehead, where the subject may be fixed, pull the cap over the head. Next, adjust the electrode on position Cz at the half way between Nasion and Inion. The frontpolar and occipital electrode should now be placed correctly. Otherwise choose another cap size. When measured correctly, the frontpolar electrodes (Fp1/Fp2) should lie directly above the eyebrows.
3. Make sure, that the position Cz is centered between the left and the right side, as well that the lateral electrode position is symmetrical.
4. Fix the cap with the chin strap and make sure that the cap fits right. The wide part of the belt should fasten at the chin, not at the neck. The belt tension should provide a reliable fixing of cap.

### Attention!

Don't pull the belt too tight. The patient should feel comfortable even the recording

## Electrodes and impedances (*MultiCap-Base*)

1. In order to minimize the impedance, the electrode position on the scalp should be prepared. We recommend the abrasive gel "EVERY" or "NUPREP" which both have very good qualities for preparation. By the user-friendly basecut of the MultiCap, all electrode positions can be reached quite simply.

2. Put the electrodes onto the desired positions. For this, press down the electrode into the electrode-body until it's be fixed.

3. Make sure, that the cables of the electrodes are not tangled or folded. Use the provided velcro-strap to focus the cables (see photo).

4. In order to raise the contact-surface, use a conductive gel like "NEURGEL", or for long time recordings the "SYNAPSE" gel that you can fill with the provided blunt needle kit into the intended gap (notch) of the electrode.

## Electrodes and impedances (*MultiCap-Soft*)

1. In order to minimize the impedance, the electrode position on the scalp should be prepared. We recommend the abrasive gel "EVERY" or "NUPREP" which both have very good qualities for preparation. By the user-friendly basecut of the MultiCap, all electrode positions can be reached quite simply.

2. Put the electrodes onto the desired positions. For this, press down the electrode into the electrode- body until it's be fixed.

3. In order to raise the contact-surface, use a conductive gel like "NEURGEL", or for long time recordings the "SYNAPSE" gel that you can fill with the provided blunt needle kit into the intended gap of the electrode.

## Electrodes and impedances (MultiCap-Cup)

1. In order to minimize the impedance, the electrode position on the scalp should be prepared. We recommend the abrasive gel "EVERY" or "NUPREP" which both have very good qualities for preparation. By the user-friendly basecut of the MultiCap, all electrode positions can be reached quite simply.

2. Put the electrodes onto the desired positions. Place the electrode on the wide part of the holder (photo)

3. In order to raise the contact-surface, use a conductive gel like "NEURGEL", or for long time recordings the "SYNAPSE" gel that you can fill with the provided blunt needle kit into the intended gap of the electrode.

## Electrodes and impedances (*MultiCap-Flat and MultiCap Doc*)

1. In order to minimize the impedance, the electrode position on the scalp should be prepared. We recommend the abrasive gel "EVERY" or "NUPREP" which both have very good qualities for preparation. By the user-friendly basecut of the MultiCap, all electrode positions can be reached quite simply.

2. This type of cap can be used without the need for additional electrode holders. Thanks to innovative technology, the electrodes are attached directly to the cap. Simple fitting makes it possible to adapt the cap quickly and individually to your needs.

3. Make sure, that the cables of the electrodes are not tangled or folded. Use the provided velcro-strap to focus the cables (see photo).

4. In order to raise the contact-surface, use a conductive gel like "NEURGEL", or for long time recordings the "SYNAPSE" gel that you can fill with the provided blunt needle kit into the intended gap of the electrode.


## Cleaning and handling of caps and electrodes

### Attention!

Electrodes and cap should be cleaned promptly after using, before the gel-resistances can be dry.

English

To remove the electrode, pull them with slow rotational motions out of the electrode-holder.

To remove the electrode, push them forward out of the holder. Pull at the cable-bend-protection, not at the cable.

To remove the electrode, pull them with slow rotational motions out of the electrode-body. Attention! Touch the electrode-body to

Stretch carefully the material of the cap to remove the electrode. Pull at the cable-bend-protection, not at the cable.

### Cleaning:

After each use, they must be cleaned in warm water and with mild detergents, e.g., Ivory, and then rinsed and dried.

Please do not use aggressive cleaning agents or shampoos.

Afterwards, rinse the hoods and the electrodes with lukewarm and clear water. Finally, dry the material in the air or with a towel.

Do NOT bring the electrodes into contact with metal.

The electrodes and hoods should be stored in a dark and dry place.

### Attention!

Improper treatment may lead to damages. Do not clean the products with abrasive cleaners or strongly disinfecting solutions which are based on, e.g., acetones.

### Disinfection:

We recommend the instrument disinfectants Gigasept FF (new) and Korsoplex plus for which the material compatibility is evidenced.

Please observe the instructions of the manufacturer of these agents.

When choosing the cleaning and disinfecting agents to be used, the material compatibility and the proven efficacy of the agents are to be considered.

(List of disinfectants – VAH)

### Attention!

Do not soak electrodes in a chloride- or saline-solution.

By this, the electrodes will corrode.

Do not use hot sterilisation methods (steam), because the wire isolation can be damaged.

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GVB-geliMED GmbH