

The needle
of the Experts!

NEW

MYOLINE
xpert premium

IMPROVED BEVEL
for pain-free insertion

NEW SURFACE FINISH
to improve the sliding characteristics of the needles

EVEN MORE FLEXIBLE CONNECTING CABLE (NEEDLE-HOLDER)
for easier handling

IMPROVED NEEDLE CONNECTION (CONTACT)
for artefact-free measurements

BETTER SIGNAL QUALITY
for high customer satisfaction

MYOLINE XP - disposable needle electrodes
Packing unit: 25 pcs.

Colour:	Length:	Diameter:	Art-No.
white	20mm	0,35mm	10009079P
violet	25mm	0,30mm	10009080P
red	30mm	0,35mm	10009081P
grey	40mm	0,35mm	10009083P
green	50mm	0,35mm	10009085P
orange	30mm	0,45mm	10009082P
black	40mm	0,45mm	10009087P
yellow	50mm	0,45mm	10009089P
blue	60mm	0,50mm	10009091P
brown	70mm	0,60mm	10009093P

CALL-UP ORDER

Use our "call-up order" option



Just let us know your estimated
annual requirement of
the product you want

We will put together a
non-binding offer for you


Test a product?

Would you like to
test a product first?

Just call us!

EMG ACCESSORIES

Compatible with all current EMG/EEG devices



GVB-gelimed E. Kfm.
Phn.: +49 45 51 - 95 67 30
Fax: +49 45 51 - 95 67 33
info@gvb-gelimed.de



www.gvb-gelimed.com



Your partner
for medical accessories

EMG ELECTRODES

Surface stimulation electrode

Dual-channel surface stimulation electrode with removable plastic holder. Each electrode has two 7mm felt stimulation points and can be secured using a Velcro strip. The connecting cable is 150cm long and has a 5-pole 270° DIN connector.

Stimulation electrode with 150cm connecting cable

and 5-channel 270° DIN connector ST90-671
and 1.5mm DIN connector ST90-626

Spare felt, 7mm round (10 pcs.) ST90-006




Shielded stimulation electrodes

Shielded stimulation electrodes with ergonomic grip for easier handling. The small electrode has two 3mm stimulation points with a gap of 10mm, the large electrode has 6mm stimulation points with a gap of 24mm. Both electrodes can also be secured with a Velcro strip.

Small electrode with 150cm cable and 5-pole 270° DIN connector and 1.5mm DIN connector SP27-671
SP27-626
Stimulation felt, 3mm round (10 pcs.) SP27-000
Stimulation balls 3 mm (1 pair) SP27-003

Large electrode with 150cm cable and 5-pole 270° DIN connector and 1.5mm DIN connector SP28-671
SP28-626
Stimulation felt, 6mm round (10 pcs.) SP28-000
Stimulation balls 6mm (1 pair) SP28-007



Shielded bar stimulation and recording electrodes

These shielded bar stimulation and recording electrodes with two 20mm x 7mm felt inserts can be used for conductance as well as stimulation. There is a choice of gap between the stimulation points of 20mm or 40mm. The connecting cable is 150 cm long and can be provided with various different connection connectors. In the basic version, 2 Velcro strips 40cm and 100cm in length are included as well as the electrode.

EMG ELECTRODES

Finger/toe stimulation and recording electrode (can be used without gel) with shielded connecting cable.

100 cm cable with 1.5 DIN connector SQ80-426
150 cm cable with 1.5 DIN connector SQ80-626
150 cm cable with 5-pole 270° connector SQ80-671
100 cm cable with 5-pole 270° connector SQ80-471



SQ50..




Electrode with standard loops made of stainless steel, with shielded connecting cable

100 cm cable with 1.5 DIN connector SQ50-426
150 cm cable with 1.5 DIN connector SQ50-626
150 cm cable with 5-pole 270° connector SQ50-471
100 cm cable with 5-pole 270° connector SQ50-671

Electrode with fine loops made of stainless steel, with shielded connecting cable

100 cm cable with 1.5 DIN connector SQ55-426
150 cm cable with 1.5 DIN connector SQ55-626
150 cm cable with 5-pole 270° connector SQ55-671
100 cm cable with 5-pole 270° connector SQ55-471



SA80..



Waffle surface electrode

Waffle recording electrode for sensitive irritation and motor response potential on extremities. Also suitable for conductance on the face. Two 8mm V2A small corrugated steel plates are securely attached to a 100cm long, highly flexible unshielded cable. White is the active electrode, brown the reference electrode.

with 5-pole DIN connector SA80-671
with 1.5 mm DIN connector SA80-626

Disc surface electrodes

10mm disc electrodes made of stainless steel, 1 pair of twisted red and black cables, 100 cm long with a choice of connectors.

with 5-pole 270° DIN connector (shielded cable) SA70-471
with 1.5 mm DIN connector SA70-426



PASTES AND GELS

A non-drying neurodiagnostic contact paste. No lumps, no second application necessary, easy to wash off, easy to clean, no residues.


Packing unit: 3 pcs. 115g jar Ten-20-SET
115g tube Ten-20-T-SET
230g jar Ten-20-8-SET



NUPREP

Modern skin preparation gel (recommended for sensitive skin). No skin irritations, pleasant for patients.

Packing unit: 3 pcs. Content: 115g NUPREP-SET



NEURGEL

Neurgel is a highly flexible and highly conductive water-soluble gel with a pH value of 6 to 7. Neurgel was developed specially for EEG, EMG and EP recordings.

250g jar GEL250V
250g bottle GEL250F



EVERI

A preparation paste specially developed to reduce skin transfer resistance. For preparation and use in EEG-EP-EMG-EKG

Content: 160ml EVERY

Specially developed adhesive conducting paste with excellent bonding properties. In combination with our Gold CUP electrode, this paste is ideal for long-term recordings.

Content: 100g 10 pcs.: SEC-2
1 pcs: EC-2

