

# CALL-UP ORDER

**GVB**  
online

The needle  
from the experts!

Use our "call-up order" option!



Just let us know your estimated  
annual requirement of the  
product you want

We will put together a  
non-binding offer for you

**Test a product?**

Would you like to test  
a product first?

Just give us a call!



Visit our online shop today:  
[www.gvb-gelimed.de](http://www.gvb-gelimed.de)

or call us on:  
+49 45 51 - 95 67 30

## Online NEWSLETTER

Register for our newsletter  
and you will always get  
the latest information  
to keep up-to-date



**MYOLINE**  
xpert premium *XP*

**NEW**

**YOUR CONTACT NO. FOR CALL-UP ORDERS  
OR TEST PRODUCTS  
+49 45 51 - 95 67 30**

GVB-gelIMED E. Kfm.  
Phn.: +49 45 51 - 95 67 30  
Fax: +49 45 51 - 95 67 33  
[info@gvb-gelimed.de](mailto:info@gvb-gelimed.de)

  
your partner for medical accessories  
[www.gvb-gelimed.com](http://www.gvb-gelimed.com)

 gelimed

your partner  
for medical accessories

### Connecting cable

100cm cable: Code  
 5-pole connector 100NHN71  
 5-pole connector (Jaeger Toennies) 100NHN78  
 1.5mm DIN connector 100NHN26

200cm cable: Code  
 5-pole connector 200NHN71  
 1.5mm DIN connector 200NHN26



### “MYOLINE XP” disposable needle electrodes

Packing unit: 25 pcs.

| Colour: | Length: | Diameter: | Art.-No.: |
|---------|---------|-----------|-----------|
| white   | 20mm    | 0,35mm    | 10009079P |
| violet  | 25mm    | 0,30mm    | 10009080P |
| red     | 30mm    | 0,35mm    | 10009081P |
| grey    | 40mm    | 0,35mm    | 10009083P |
| green   | 50mm    | 0,35mm    | 10009085P |
| orange  | 30mm    | 0,45mm    | 10009082P |
| black   | 40mm    | 0,45mm    | 10009087P |
| yellow  | 50mm    | 0,45mm    | 10009089P |
| blue    | 60mm    | 0,50mm    | 10009091P |
| brown   | 70mm    | 0,60mm    | 10009093P |

**IMPROVED BEVEL**  
for pain-free insertion

**NEW SURFACE QUALITY**  
to improve the sliding characteristics of the needles

**EVEN MORE FLEXIBLE CONNECTING CABLE (NEEDLE-HOLDER)**  
for easier handling

**IMPROVED NEEDLE CONNECTION (CONTACT)**  
for artefact-free measurements

**BETTER SIGNAL QUALITY**  
for high customer satisfaction



The needle  
of the **Xperts!**