

NEW PRODUCTS!

SD57



NEW

Disposable subdermal steel needle with concentric bevel

13 mm needle length
0,20mm needle diameter.

SD57-426-M

Set of 25: 5 x 5 needles in 5 colours, 100cm cable, 1,5mm DIN connector

SD51-007-E



Detail:
Needle connection

NEW



Connection cable 10020671
(3-channel also available)

SD51-007-E

With 0.7mm PIN connection, for connecting to cable 10020671. The cable (length 100cm or 200cm) can be fitted with either a 5-pole DIN connector or a 1.5mm DIN connector

CALL-UP ORDER for subdermal needles

Use our "call-up order" option!



Just let us know your estimated annual requirement of the product you want.

We will put together a non-binding offer for you

Test the product?

Would you like to test a product first?

Just give us a call!

EBG EMG EP IOM
NEURO.DART.
SUBDERMAL NEEDLE ELECTRODES

SUBDERMAL NEEDLES



You can also find all needle electrodes and other accessories online at www.gvb-gelimed.com

GVB-gelimed E. Kfm.
Phn.: +49 45 51 - 95 67 30
Fax: +49 45 51 - 95 67 33
info@gvb-gelimed.de


your partner for medical accessories
www.gvb-gelimed.com



your partner for
medical accessories

EEG EMG EP ICM
NEURO.DART.
SUBDERMAL NEEDLE ELECTRODES

UNIMPEDED RECORDINGS SUITED TO STIMULATION PARTICULARLY SHARP BEVEL

SD52-426-M



Subdermal needle with 100cm twisted connecting cable and 1.5mm DIN connector
13mm needle length / 0.4mm needle diameter

Packing unit: 25 pairs in 5 colours

SD51...



Single-use Subdermal steel needles with concentric bevel. 13mm needle length / 0.40mm needle diameter. Options: 100cm, 150cm or 250cm connecting cables.

Packing unit: Set of 5 x 5 needles in 5 colours

SD51M...



Single-use subdermal steel needles with concentric bevel and soft rubber needle connection. 13mm needle length / 0.40mm needle diameter. Options: 100cm, 150cm or 250cm connecting cables.

Packing unit: Set of 5 x 5 needles in 5 colours

OTHER SUBDERMAL NEEDLES IN OUR RANGE:

SD55...



Single-use corkscrew electrode for recording EEG and SEP signals.
A choice of 100cm, 150cm or 250cm cables and 1.5mm DIN connector.

Packing unit: 25 pcs.

MN3512F110-S



Single-use one-channel subdermal steel needles with 1mm connection, needle length 12mm, diameter 0.35mm. A choice of 100cm or 150cm connecting cable.

Packing unit: 5 x 5 pcs.

SD51-111



Single-use subdermal steel needles, needle length 12mm, 0.35mm diameter, 1mm PIN connector, 8cm cable.

Packing unit: 5 x 5 pcs.

MN3512H150



Single-use subdermal hook needle with 8cm cable and 1mm PIN connector.

Packing unit: 5 pcs.