

LUER DISPOSABLE FILTER



Model 14809

This low cost hydrophobic disposable filter allow you to use any type pf cannula with the pressure sensor, eliminating the fear of contaminating the sensor and patients cross contamination.

SPECIFICATIONS

Filter specifications	Cannula connection:
	Female 4mm Luer
	Sensor connection:
	Male 4mm Luer.
	Single use only.
Physical Size	1" (25m) outside diameter
Weight	0.2 oz. (5 g)

ORDERING INFORMATION

14809	Single use Luer filter.
-------	-------------------------

SleepSense®
Sensors that make sense